



Belle® can be located with **WiFi** location technologies.

It is **lightweight** and easy to wear or carry with the included lanyard and belt clip.

Belle® lasts up to **30 days** per charge.

It comes with a **simple charging cradle** to make recharging easy.



In addition to a lanyard, Belle also comes with a **snap-on belt clip** for maximum flexibility as you wear or carry the pendant.

Belle is **IP67 water-resistant**, though it should not be submerged.




Order Belle today and enjoy the freedom to explore and the confidence to do what you love.

4G LTE

belle®

go with confidence!



mobile
medical alert
system



Belle® is a small, easy to use mobile medical alert pendant that offers you confidence to live life on your terms.

It is certified on the Verizon 4G LTE network, so you can easily get help at home or away from home with the press of a button, across the US.

Here's how it works...

speaker

With this powerful speaker, it's easy to hear the specialist answer your call.

call button

Just press one button to speak with a specialist and get help. It's that simple.

circular light

The light around the button lets you know when a call is in progress and when the pendant needs to be charged.



microphone

A powerful microphone allows you to speak with a kind, courteous specialist. They can send emergency help or contact a loved one to assist.