



Belle® is a registered trademark of Freeus, LLC.

4G LTE

belle®

go with confidence!

mobile  
medical alert  
system

*Belle® is a small, easy to use mobile medical alert pendant that offers you confidence to live life on your terms. It is certified on the Verizon 4G LTE network, so you can easily get help at home or away from home with the press of a button, across the US.*

## Here's how it works...

### speaker

With this powerful speaker, it's easy to hear the specialist answer your call.

### call button

Just press one button to speak with a specialist and get help. It's that simple.

### circular light

The light around the button lets you know when a call is in progress and when the pendant needs to be charged.

### battery status button

Press anytime for a real-time battery level update

### microphone

A powerful microphone allows you to speak with a kind, courteous specialist. They can send emergency help or contact a loved one to assist.



Belle® lasts up to **30 days** per charge, and it comes with a **simple charging cradle** to make recharging easy. Belle is lightweight and easy to wear or carry with the **included lanyard and belt clip**.

It includes **WiFi location technologies**. It is IP67 **water-resistant** but should not be submerged.



**HelpLine®**

Medical Alert Systems

207-892-1454

Home Support Services, LLC  
744 Roosevelt Trail, Suite 210  
Windham, Maine 04062  
[www.helplinemedicalalert.com](http://www.helplinemedicalalert.com)